



LightCOM

Lighthouse Consortium on
Obesity Management

Study name: **Lighthouse Consortium on Obesity Management**

Short name: **LightCOM**

SOP title: **Standard Operating Procedure for Patient Care-DEXA Scanner**

SOP short title: **SOP on DXA**

SOP Version: **1**

Date: **09.04.2024**

Authored by:

Odette Dawkins

Date:

25.09.2023

TABEL OF CONTENTS

1	PURPOSE	3
2	INTRODUCTION	3
3	SCOPE	3
4	RESPONSIBILITIES	3
4.1	Chief Investigator/Principle Investigator or delegate	3
5	SPECIFIC PROCEDURE	3
5.1	Equipment Check.....	3
6	FORMS/TEMPLATES TO BE USED	4
7	INTERNAL AND EXTERNAL REFERENCES	4
7.1	Internal References.....	4
7.1	External References	4
8	HISTORICAL RECORD OF SOP	5

1 PURPOSE

The purpose of this Standard Operating Procedure (SOP) is to describe the standard procedures to be followed when asking a participant to undertake a DEXA Scan during the LightCOM studies which includes LightBAR, LightWAY and LightCARE.

2 INTRODUCTION

A DEXA Scan is a Dual Energy X-ray Absorptiometry to test the bone density of participants. It will be conducted on participants at two timepoints within the study. The scan results will be held on the study database.

3 SCOPE

This SOP applies to clinical research carried out under the LightCOM research studies which incorporates the individual studies of LightBAR, LightWAY and LightCARE. The sponsor is Department of Medicine, Copenhagen University Hospital Amager and Hvidovre, Copenhagen, Denmark.

4 RESPONSIBILITIES

4.1 Chief Investigator/Principle Investigator or delegate

Will be responsible for ensuring the staff member assigned to use the DEXA Scanner is suitably qualified to undertake the procedure and has received adequate and up to date training.

5 SPECIFIC PROCEDURE

5.1 Equipment Check

- Test the emergency 'Stop' button on the scanner and ensure it is in full working order.
- Make sure the X-ray and Shutter lights are off and that a message on the computer monitor indicates the emergency stop button is activated.
- Ensure there is a sign displayed on the outer door of the scanner room and is in full working order to show when scanning is in progress. Check the outer door can be securely closed.
- Follow the local quality assurance protocols to calibrate the machine prior to commencing any participant scans.
- Clean the scanning table in line with local infection control protocols.

5.2 Patient Care

- UK: Check the participant has received the latest version of the local Patient Information Sheet for a DEXA Scan prior to the scan date and is fully aware of the scan procedure.
- Ensure full signed consent has been obtained and complete the scan log in line with local procedures.

- If the participant to be scanned is female and under the age of 55yrs, ask if they think they may be pregnant. If unsure, ask permission to conduct a pregnancy test using a urine sample test and record the results in the scan log.
- Ask the participant to empty their bladder
- Ask the participant to remove any items of jewellery or clothing which contain metal and if required, provide them with a hospital gown. Ask the participant to remove their shoes.
- Check the participant's DOB, Study ID and record their height and weight in the clinical system.
- Ask the participant to lie on the scanning table and remain still during the procedure.

5.3 Scan Procedure

- Follow the local specified scan instructions during the scan procedure. Select 'Total Body' measure and check the position of the participant on the scanning couch, adjust if necessary.
- Start the scan following the specified guidelines for the scan equipment. Do not enter the scanning area until the scan is over.
- The full body scan should last approximately 15 minutes. Regional bone scan of hip, lower back, and non-dominant forearm (forearm measurements for LightBAR only) last approximately 30-40 min, however, timings may differ by machine.
- For hip and forearm, it is important to register the side scanned (L vs R) in the eCRF and the same side used for the follow up measure.
- Once the procedure is completed, ask the participant to carefully leave the scanning table ensuring they are safe to do so after being in the supine position.
- Enter the scan results onto the study database as per Study Protocol.

6 FORMS/TEMPLATES TO BE USED

Use Study database in OpenClinica for DEXA Scan results as per study protocol.

7 INTERNAL AND EXTERNAL REFERENCES

7.1 Internal References

Design and Development of SOP – CTRG Template SOP V3.2

7.2 External References

NHS OBB DEXA Scan

Surry and Sussex NHS Trust DEXA Scan

8 HISTORICAL RECORD OF SOP

Rev# / Date	Written/Revised By	Revision Description	Approval
0. Written	Odette Dawkins		
1. Edited			
2.			
3.			
4.			
5.			
6.			
Removed / replaced:			